

Principle Behind the CO₂ Correction Factor

The BreathTrackers SC and H2+ utilize a CO₂ correction factor technique to minimize error caused by improper sampling techniques. The CO₂ correction factor is based on the concept that carbon dioxide is present in alveolar (lung) air at a virtually constant concentration; while CO₂ in room air is virtually zero (in fact it is present, but in extremely trace concentrations). Therefore, if an alveolar air sample is accidentally contaminated (mixed) with room air, the CO₂ concentration in the sample will be reduced, as will any other trace gases in the sample (in our case H₂ and/or CH₄). By knowing the degree to which the CO₂ is diluted, it is possible to apply a correction to the analysis of each trace gas of interest (again in our case H₂ and/or CH₄), and be able to estimate the “true alveolar” concentration of these trace gases. The sample concentrations of H₂ and/or CH₄ are multiplied by the factor calculated from:

$$\text{FACTOR} = \text{Alveolar CO}_2 \text{ concentration} / \text{Sample CO}_2 \text{ concentration}$$

CO₂ is the physiological regulator of breathing, and the whole breathing system is dedicated to keeping the alveolar CO₂-pressure (PACO₂) constant at 40 mm Hg (torr). Therefore CO₂ is the most reliable “normalizing” component in the sample, because it ordinarily has the most constant alveolar composition of any gas in the sample.

Alveolar PCO₂ remains constant at 40 torr among normal individuals if ventilation is normal. The percent of CO₂ in an alveolar sample is affected by the barometric pressure (altitude) at which the sample is collected. Alveolar air with a PCO₂ of 40 torr in Miami (at sea level) will have a CO₂ concentration of about 5.5% in dry air (40/(760-37)), while alveolar air in Denver (where barometric pressure is closer to 625 torr) will be near 6.8% (40/(625-37)). Significant differences in barometric pressure exist at different altitudes, as demonstrated by Miami and Denver. However, using a single correction factor, alveolar concentration will simplify the process without introducing significant error, because all the samples will be normalized to the same (constant) CO₂ level.

Research studies, using an alveolar concentration of 5.5% will be adequate for calculating the correction factor for CO₂. However when the absolute alveolar pressure for the trace gases is important, you may adjust the instrument to your specification.

CO₂ (Correction Feature) for Methane (CH₄) is only available on the BreathTracker SC.

