

# Collection System Component Matrix

THE COLLECTION SYSTEM COMPONENT MATRIX CAN BE USED TO HELP END-USERS DETERMINE THE PRODUCT COMPONENTS UTILIZED IN EACH COLLECTION SYSTEM.

## COLLECTION SYSTEMS

CATALOG # & DESCRIPTION	GaSampler™	KidSampler™	BabySampler™	Additional Info.
QT00830-P - 750mL Single-Patient Collection Bag	X	X		<b>Includes:</b> Large port cap and one-way stopcock
QT00834-P - 250mL Single-Patient Collection Bag	X	X	X	<b>Includes:</b> Large port cap and one-way stopcock
QT00841-P - 720mL Multi-Patient Collection Bag	X	X		<b>Includes:</b> Large port cap and one-way stopcock
QT00844-P - 250ml Multi-Patient Collection Bag	X	X	X	<b>Includes:</b> Large port cap and one-way stopcock
QT00842-P - 250mL Multi-Patient Sample Holding Bag	X	X	X	<b>Includes:</b> One-way stopcock and male-male luer adapter
QT00843-P - 400mL Discard Bag	X	X		
QT00850-P - Tee-Piece Only	X	X		
QT00851-P - Flutter Valve for Tee-Piece	X	X		
QT00854-P - Tee-Mouthpiece Assembly	X	X		<u>Product Assembled With:</u> 1-QT00850-P, 1-QT00851-P, 1-QT00991-P
QT00855-P - Disposable Mask Adapter for Pediatric Face Masks			X	Only needed for Infant and Neonatal Face Masks, attaches face mask to Tee-Piece.
QT00859-P - Low Deadspace - One-Way Tee-Piece Assembly			X	
QT00881-L - 5.25cm, Neonatal Face Mask			X	
QT00882-L - 6.0cm, Infant Face Mask			X	
QT00883-L - 7.0cm, Toddler Face Mask		X	X	
QT00884-L - 8.5cm, Child Face Mask	X	X	X	
QT00885-L - 9.0cm, Adult Face Mask	X	X		
QT00991-P - Universal Plastic Mouthpiece	X	X		
QT01727-V - One-Way Plastic Stopcock	X	X	X	

Only products sold as “Multi-Patient” can be safely used inter-patient pending proper use, storage, handling and cleaning practices.

*You should follow your institution's practices regarding disposable/inter-patient products.*

All other products not listed as Multi-Patient are for Single-Patient Use and are not recommended for reuse with different patients.